

LOT: 12302-52
EXP: 2023-07-23
Open Vial Stability: 16 days

EXIGO		Mean	± Range	Expected Range	
WBC	10 ⁹ /L	8.7	0.6	8.1	9.3
Hct	%	30.2	2.5	27.7	32.7
MCV	fL	73.2	5.0	68.2	78.2
RBC	10 ¹² /L	4.13	0.18	3.95	4.31
Hgb	g/dL	10.2	0.4	9.8	10.6
MCH	pg	24.7	1.5	23.2	26.2
MCHC	g/dL	33.8	2.2	31.6	36.0
PLT	10 ⁹ /L	213	30	183	243

This assay sheet is also available on :
www.exigo-vet.com

For EXIGO/EXIGO EOS with barcode reader

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Instr.: Exigo
Lot: 12302-52
Prod.: Vet-Con Norm
ID: 1230252+
Exp.: 23-Jul-2023

