

LOT: 12111-52
EXP: 2022-04-24
Open Vial Stability: 16 days

| Medonic CA530, CA620 | | Mean | ± Range | Expected Range | |
|----------------------|---------------------|------|---------|----------------|------|
| WBC | 10 ⁹ /L | 8.6 | 0.6 | 8.0 | 9.2 |
| Hct | % | 33.7 | 2.5 | 31.2 | 36.2 |
| MCV | fL | 80.7 | 5.0 | 75.7 | 85.7 |
| RBC | 10 ¹² /L | 4.18 | 0.18 | 4.00 | 4.36 |
| Hgb | g/dL | 11.5 | 0.4 | 11.1 | 11.9 |
| MCH | pg | 27.5 | 1.5 | 26.0 | 29.0 |
| MCHC | g/dL | 34.1 | 2.2 | 31.9 | 36.3 |
| PLT | 10 ⁹ /L | 232 | 30 | 202 | 262 |

| EXIGO | | Mean | ± Range | Expected Range | |
|-------|---------------------|------|---------|----------------|------|
| WBC | 10 ⁹ /L | 8.3 | 0.6 | 7.7 | 8.9 |
| Hct | % | 34.5 | 2.5 | 32.0 | 37.0 |
| MCV | fL | 83.2 | 5.0 | 78.2 | 88.2 |
| RBC | 10 ¹² /L | 4.15 | 0.18 | 3.97 | 4.33 |
| Hgb | g/dL | 11.6 | 0.4 | 11.2 | 12.0 |
| MCH | pg | 28.0 | 1.5 | 26.5 | 29.5 |
| MCHC | g/dL | 33.6 | 2.2 | 31.4 | 35.8 |
| PLT | 10 ⁹ /L | 207 | 30 | 177 | 237 |

This assay sheet is also available on :

www.medonic.se
www.exigo-vet.com

For EXIGO/EXIGO EOS with barcode reader

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Instr.: Exigo
Lot: 12111-52
Prod.: Vet-Con Norm
ID: 1211152+
Exp.: 24-Apr-2022

