

LOT: 12211-52
EXP: 2023-04-23
Open Vial Stability: 16 days

EXIGO		Mean	± Range	Expected Range	
WBC	10 ⁹ /L	8.6	0.6	8.0	9.2
Hct	%	31.3	2.5	28.8	33.8
MCV	fL	75.4	5.0	70.4	80.4
RBC	10 ¹² /L	4.15	0.18	3.97	4.33
Hgb	g/dL	10.6	0.4	10.2	11.0
MCH	pg	25.5	1.5	24.0	27.0
MCHC	g/dL	33.9	2.2	31.7	36.1
PLT	10 ⁹ /L	220	30	190	250

This assay sheet is also available on :
www.exigo-vet.com

For EXIGO/EXIGO EOS with barcode reader

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Instr.: Exigo
Lot: 12211-52
Prod.: Vet-Con Norm
ID: 1221152+
Exp.: 23-Apr-2023

