

LOT: 12007-52
EXP: 2020-12-21
Open Vial Stability: 16 days

Medonic CA530, CA620		Mean	± Range	Expected Range	
WBC	10 ⁹ /L	8,6	0,6	8,0	9,2
Hct	%	32,3	2,5	29,8	34,8
MCV	fL	78,8	5,0	73,8	83,8
RBC	10 ¹² /L	4,10	0,18	3,92	4,28
Hgb	g/dL	11,2	0,4	10,8	11,6
MCH	pg	27,3	1,5	25,8	28,8
MCHC	g/dL	34,7	2,2	32,5	36,9
PLT	10 ⁹ /L	242	30	212	272

EXIGO		Mean	± Range	Expected Range	
WBC	10 ⁹ /L	7,8	0,6	7,2	8,4
Hct	%	32,3	2,5	29,8	34,8
MCV	fL	81,1	5,0	76,1	86,1
RBC	10 ¹² /L	3,98	0,18	3,80	4,16
Hgb	g/dL	11,3	0,4	10,9	11,7
MCH	pg	28,4	1,5	26,9	29,9
MCHC	g/dL	35,0	2,2	32,8	37,2
PLT	10 ⁹ /L	217	30	187	247

This assay sheet is also available on :

www.medonic.se
www.exigo-vet.com

For EXIGO/EXIGO EOS with barcode reader

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Instr.: Exigo
Lot: 12007-52
Prod.: Vet-Con Norm
ID: 1200752+
Exp.: 21-Dec-2020

