

# Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days <sup>1</sup>



Control	Lot	
Tri-Level	22201-1K	2022-06-07
Low	22201-11	2022-06-07
Normal	22201-12	2022-06-07
High	22201-13	2022-06-10

MEDONIC		CONTROL	L
CA620/530		LOT	22201-11
		Mean <sup>2</sup>	± Range <sup>3</sup>
RBC	10 <sup>12</sup> /L	2.27	0.12
MCV	fL	63.5	5.0
RDW	%	14.6	4.0
Hct	%	14.4	2.5
Plt	10 <sup>9</sup> /L	88	15
MPV	fL	9.8	1.5
WBC	10 <sup>9</sup> /L	3.5	0.3
Hgb	g/dL	5.6	0.3
MCH	pg	24.7	1.5
MCHC	g/dL	38.9	3.0
Lymph#	10 <sup>9</sup> /L	1.6	0.4
Mid#	10 <sup>9</sup> /L	0.2	0.3
Gran#	10 <sup>9</sup> /L	1.7	0.8
Lymph	%	45.2	7.0
Mid	%	5.2	5.0
Gran	%	49.6	7.0

CONTROL		N	
LOT		22201-12	
		Mean <sup>2</sup>	± Range <sup>3</sup>
		4.13	0.18
		76.8	5.0
		14.5	4.0
		31.7	3.0
		249	30
		10.4	1.5
		8.9	0.6
		11.5	0.4
		27.8	1.5
		36.3	3.0
		4.3	0.6
		0.6	0.5
		4.0	1.0
		48.4	5.0
		6.7	5.0
		44.9	5.0

CONTROL		H	
LOT		22201-13	
		Mean <sup>2</sup>	± Range <sup>3</sup>
		5.00	0.22
		89.6	5.0
		14.2	4.0
		44.8	3.5
		555	60
		11.5	1.5
		21.5	1.8
		15.9	0.5
		31.8	1.5
		35.5	3.0
		13.3	1.5
		1.3	1.0
		6.9	2.5
		61.9	5.0
		6.1	5.0
		32.0	5.0

Medonic M-series/M32		Mean <sup>2</sup>	± Range <sup>3</sup>
RBC	10 <sup>12</sup> /L	2.26	0.12
MCV	fL	65.0	5.0
RDW	%	14.7	4.0
Hct	%	14.7	2.5
Plt	10 <sup>9</sup> /L	79	15
MPV	fL	10.0	1.5
WBC	10 <sup>9</sup> /L	3.5	0.3
Hgb	g/dL	5.5	0.3
MCH	pg	24.3	1.5
MCHC	g/dL	37.4	3.0
Lymph#	10 <sup>9</sup> /L	1.5	0.4
Mid#	10 <sup>9</sup> /L	0.2	0.3
Gran#	10 <sup>9</sup> /L	1.8	0.8
Lymph	%	43.9	7.0
Mid	%	5.0	5.0
Gran	%	51.1	7.0

	Mean <sup>2</sup>	± Range <sup>3</sup>
	4.14	0.18
	79.0	5.0
	14.4	4.0
	32.7	3.0
	228	30
	10.4	1.5
	8.5	0.6
	11.6	0.4
	28.0	1.5
	35.5	3.0
	4.1	0.6
	0.5	0.5
	3.9	1.0
	47.7	5.0
	6.3	5.0
	46.0	5.0

	Mean <sup>2</sup>	± Range <sup>3</sup>
	5.06	0.22
	93.0	5.0
	13.8	4.0
	47.1	3.5
	509	60
	11.3	1.5
	20.6	1.8
	16.2	0.5
	32.0	1.5
	34.4	3.0
	12.6	1.5
	1.2	1.0
	6.8	2.5
	61.4	5.0
	5.7	5.0
	32.9	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
  - hr Stabilnost otvorene bočice 14 dani (hr)
  - cs Stabilita otevřené ampulky 14 dní (cs)
  - da Åben rør stabilitet 14 dage (da)
  - et Stabiilsus avatud vialis 14 päevade (et)
  - fr Stabilité en flacon ouvert 14 jour (fr)
  - de Stabilität geöffneter Flaschen 14 tage (de)
  - el Διατήρηση δέγματος μετά την αποσφράγιση 14 ημέρα (el)
  - it Stabilità della fiala aperta 14 giorni (it)
  - lv Atvērt flakonu stabilitāti 14 dienas (lv)
  - lt Stabilumas atidarius buteliuką 14 dienas (lt)
  - no Åpen rør/glass stabilitet 14 dager (no)
  - pl Trwałość po otwarciu fiołki 14 Liczba dni (pl)
  - pt Estabilidade após abertura do frasco 14 dias (pt)
  - ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
  - ru Стабильность вскрытого флакона 14 дней (ru)
  - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
  - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
  - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
  - es Estabilidad de la cápsula abierta 14 días (es)
  - sv Hållbarhet för öppnad flaska 14 dagar (sv)
  - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Sredno
  - Sredstvo
  - Průměr
  - Middelverdi
  - Keskmine
  - Moyenne
  - Mittelwert
  - Mέση
  - Media
  - Nozīmē
  - Vidurkis
  - Middelverdi
  - Średnia
  - Média
  - Medie
  - Среднее
  - Srednja vrednost
  - Priemer
  - Povprečno
  - Media
  - Medelvärde
  - Ortalama

- 3
- Обхват
  - Raspon
  - Rozsah
  - Område
  - Vahemik
  - Intervalle
  - Bereich
  - Αναμενόμενο εύρος α
  - Intervallo
  - Diapazons
  - Intervalas
  - Område
  - Zakres
  - Intervalo
  - Interval
  - диапазон
  - Opseg
  - Rozsah
  - Območje
  - Intervalo
  - Intervall
  - Aralığı

25121-4

1504054



[http://www.medonic.se/  
support/downloads](http://www.medonic.se/support/downloads)

  
Boule Medical AB  
Domnarvsgatan 4  
SE-163 53 Spånga, Sweden

# For Medonic M-series / M32 with barcode reader

## M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

## M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series  
 Lot: 22201-11  
 Prod.: Con-Diff Low  
 ID: 2220111+  
 Exp.: 7-Jun-2022

Instr.: M-series  
 Lot: 22201-12  
 Prod.: Con-Diff Norm  
 ID: 2220112+  
 Exp.: 7-Jun-2022

Instr.: M-series  
 Lot: 22201-13  
 Prod.: Con-Diff High  
 ID: 2220113+  
 Exp.: 10-Jun-2022

