

Boule Con-Diff

CONTROL

IVD



Control	Lot	
Tri-Level	22301-1K	2023-06-01
Low	22301-11	2023-06-01
Normal	22301-12	2023-06-06
High	22301-13	2023-06-08

Open vial stability 14 days ¹

MEDONIC		CONTROL	L
Medonic M-series/M32		LOT	22301-11
		Mean ²	± Range ³
RBC	10 ¹² /L	2.21	0.12
MCV	fL	62.6	5.0
RDW	%	14.4	4.0
Hct	%	13.8	2.5
Plt	10 ⁹ /L	80	15
MPV	fL	10.0	1.5
WBC	10 ⁹ /L	3.3	0.3
Hgb	g/dL	5.3	0.3
MCH	pg	24.0	1.5
MCHC	g/dL	38.4	3.0
Lymph#	10 ⁹ /L	1.4	0.4
Mid#	10 ⁹ /L	0.2	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	43.9	7.0
Mid	%	5.6	5.0
Gran	%	50.5	7.0

CONTROL	N	
LOT	22301-12	
	Mean ²	± Range ³
	4.13	0.18
	83.4	5.0
	13.5	4.0
	34.4	3.0
	233	30
	10.2	1.5
	8.4	0.6
	12.3	0.4
	29.8	1.5
	35.8	3.0
	3.9	0.6
	0.5	0.5
	4.0	1.0
	46.2	5.0
	6.3	5.0
	47.5	5.0

CONTROL	H	
LOT	22301-13	
	Mean ²	± Range ³
	5.03	0.22
	90.0	5.0
	14.2	4.0
	45.3	3.5
	527	60
	10.4	1.5
	20.9	1.8
	15.8	0.5
	31.4	1.5
	34.9	3.0
	12.7	1.5
	1.2	1.0
	7.0	2.5
	60.8	5.0
	5.6	5.0
	33.6	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
 - hr Stabilnost otvorene bočice 14 dani (hr)
 - cs Stabilita otevřené ampulky 14 dní (cs)
 - da Åben rør stabilitet 14 dage (da)
 - et Stabiilsus avatud vialis 14 päevade (et)
 - fr Stabilité en flacon ouvert 14 jour (fr)
 - de Stabilität geöffneter Flaschen 14 tage (de)
 - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
 - it Stabilità della fiala aperta 14 giorni (it)
 - lv Atvērt flakonu stabilitāti 14 dienas (lv)
 - lt Stabilumas atidarius buteliuką 14 dienas (lt)
 - no Åpen rør/glass stabilitet 14 dager (no)
 - pl Trwałość po otwarciu fiołki 14 Liczba dni (pl)
 - pt Estabilidade após abertura do frasco 14 dias (pt)
 - ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
 - ru Стабильность вскрытого флакона 14 дней (ru)
 - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
 - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
 - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
 - es Estabilidad de la cápsula abierta 14 días (es)
 - sv Hållbarhet för öppnad flaska 14 dagar (sv)
 - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Средно
 - Sredstvo
 - Průměr
 - Middelværdi
 - Keskmine
 - Moyenne
 - Mittelwert
 - Μέση
 - Media
 - Nozīmē
 - Vidurkis
 - Middelverdi
 - Średnia
 - Média
 - Medie
 - Среднее
 - Srednja vrednost
 - Priemer
 - Povprečno
 - Media
 - Medelvärde
 - Ortalama
- 3
- Обхват
 - Raspon
 - Rozsah
 - Område
 - Vahemik
 - Intervalle
 - Bereich
 - Αναμενόμενο εύρος α
 - Intervallo
 - Diapazons
 - Intervalas
 - Område
 - Zakres
 - Intervalo
 - Interval
 - диапазон
 - Opseg
 - Rozsah
 - Območje
 - Intervall
 - Aralığı

<http://www.medonic.se/support/downloads>

 **Boule**
Boule Medical AB
Domnarvsgatan 4
SE-163 53 Spånga, Sweden

25121-5

1504054

For Medonic M-series / M32 with barcode reader

M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series
Lot: 22301-11
Prod.: Con-Diff Low
ID: 2230111+
Exp.: 1-Jun-2023

Instr.: M-series
Lot: 22301-12
Prod.: Con-Diff Norm
ID: 2230112+
Exp.: 6-Jun-2023

Instr.: M-series
Lot: 22301-13
Prod.: Con-Diff High
ID: 2230113+
Exp.: 8-Jun-2023

